SU KALİTE KONTROL LABORATUVARI

İSTANBUL BÜYÜKŞEHİR BELEDİYESİ

İ S K İ

İSTANBUL SU VE KANALİZASYON İDARESİ

İÇİNDEKİLER

1. BÖLÜM

SU VE KALİTENİN ADRESİ: İSKİ

2. BÖLÜM

MUSLUKLARDAN AKAN SUYUN SERÜVENİ

3. BÖLÜM

İÇME SUYUNUN NİTELİKLERİ

4. BÖLÜM

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©90 312 418 82 00 - www.turkak.org.tr

AKREDİTE OLAN PARAMETRELER

FİZİKSEL ANALİZ (Physica Analysis)		Metod ve Standart Adı	TS 226 Nisan-05
pH		SM (2005) 4500 - H+ .B. Electrometric Method	6,5-9,5
Bulanıklık (Turbidity)	NTU	SM (2005) Turbidity 2130 B. Nephelometric Method	1
KİMYASAL ANALİZ (CHEMICAL ANALYSIS)			
Amonyak (Ammonia)	mg/L	S.M 4500-NH3-D Ammonia-Selective Elektrode Method	0,5
Amonyum (Ammonium)	mg/L	TS EN ISO 14911 Nisan 2000 - Ion Chromatography	0,5
Bromat (Bromate)	mg/L	SM (2005) 4110 B-Ion Chromatography with Chemical Suppression of Eluant Conductivity	0,01
Bromür (Bromide)	mg/L	SM (2005) 4110 B-Ion Chromatography with Chemical Suppression of Eluant Conductivity	
Florür (Fluoride)	mg/L	SM (2005) 4110 B-Ion Chromatography with Chemical Suppression of Eluant Conductivity	1,5
Fosfat (Phosphate)	mg/L	SM (2005) 4110 B-Ion Chromatography with Chemical Suppression of Eluant Conductivity	
Kalsiyum (Calcium)	mg/L	TS EN ISO 14911 Nisan 2000-Ion Chromatography	
Klorür (Chloride)	mg/L	SM (2005) 4110 B-Ion Chromatography with Chemical Suppression of Eluant Conductivity	250
Lityum (Lithium)	mg/L	TS EN ISO 14911 Nisan 2000-Ion Chromatography	
Magnezyum (Magnesium)	mg/L	TS EN ISO 14911 Nisan 2000 - Ion Chromatography	
Nitrat (Nitrate)	mg/L	SM (2005) 4110 B-Ion Chromatography with Chemical Suppression of Eluant Conductivity	50
Nitrit (Nitrite)	mg/L	SM (2005) 4110 B-Ion Chromatography with Chemical Suppression of Eluant Conductivity	0,5
Potasyum (Potassium)	mg/L	TS EN ISO 14911 Nisan 2000-Ion Chromatography	
Sodyum (Sodium)	mg/L	TS EN ISO 14911 Nisan 2000-Ion Chromatography	200
Sülfat (Sulfate)	mg/L	SM (2005) 4110 B-Ion Chromatography with Chemical Suppression of Eluant Conductivity	250
Serbest Klor (Free Chlorine)	mg/L	SM (2005) 4500-Cl F DPD Ferrous Titrimetric Method	
BAKTERİYOLOJİK ANALİZ (BACTERIOLOGICAL ANALYSIS)			
T.Koli.Bakt. (T.Coli.Bact.)	100mL'de	SM (2005) 9222 B Total Coliform Membrane Filter Procedure	0
E.coli (E coli)	100mL'de	SM (2005) 9222 B Total Coliform Membrane Filter Procedure	0
T.Koli.Bakt. (T.Coli.Bact.)	EMS/100mL	SM(2005) 9221 B-C-E Total Coliform Fermentation Technique	0
Fekal Koliform (Fecal Coliform)	EMS/100mL	SM(2005) 9221 B-C-E Fecal Coliform Procedure	0
Toplam Bakteri Sayısı (HPC) 35 °C'de	kob/mL	SM(2005) 9215 A-B Pour Plate Method	-
Toplam Bakteri Sayısı (HPC) 22 °C'de	kob/mL	SM(2005) 9215 A-B Pour Plate Method	-

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


İSTANBUL BÜYÜKŞEHİR BELEDİYESİ
İ S K İ
İSTANBUL SU VE KANALİZASYON İDARESİ

SERVICE UNITS

- CHEMISTRY LABORATORY TECHNICAL CHIEFTAINCY
- BACTERIOLOGY LABORATORY TECHNICAL CHIEFTAINCY
- DISINFECTION AND DISCHARGE TECHNICAL CHIEFTAINCY
- ADMINISTRATIVE AFFAIRS DEPARTMENT CHIEFTAINCY
- SAMPLE REGISTRATION UNIT

TRAINING SERVICES

- TS EN ISO/IEC 17025/2005
"General Conditions for the Competency of
Test and Calibration laboratories" Training
- Method Validation Training
Method Validation in Bacteriological Experiments
Method Validation in Chemical Experiments
- Metrology and Calibration of Glass Materials
Training
- Inter-institution Management System Inspector
training

CONTACT ADDRESS OF DIRECTORY OF WATER QUALITY CONTROL

- Water Quality Control Director: Şahin ÖZAYDIN
Phone: 0212 231 80 28
e-mail : sozaydin@tski.gov.tr
Tel: 0212 233 23 63 / 234 50 73
- 5605 - Administration (Price Information and Notification
of Complaints)
- 5610 - To receive the results of Chemical Analysis
- 5621 - Sample Registration
- 5617 - To receive the results of Bacteriological Analysis
Fax: 0212 233 07 04
Address: Paşa mh. Yay meydanı cd. No: 2
Feriköy/Şişli - Pk. 34377

OUR QUALITY VISION

İSKİ pioneers in the development of water quality control laboratories, in accordance with national sustainable development targets and parallel to the improvements in the world.

In this context, by performing research and application projects directed to middle- and long-term requirements of our country within the scope of water quality, contributing to the economy of our country as well as contributing to science and technology.

İSKİ Water Quality Control Laboratory has aimed at providing service of good quality to its customers that meets the standard requirements of EN ISO/IEC 17025 Test and General Requirements for the Adequacy of Calibration Laboratories.

- Providing a laboratory environment that meets the requirements of EN ISO/IEC 17025 International Standard, delivering a service to its customers with ever increasing quality through a good professional and technical application,


- Becoming a nationally and internationally leading laboratory regarding reliability, honesty, impartiality, secrecy, continuous improvement and being scientific with the mission that it undertakes and with the quality of the services that it produces,

- Aiming to become a leading and competent research and control centre, performing researches in the field of Water Quality Control.

Our most important values in our institutional success are our service quality, human resources and experiences.

ACCREDITED PARAMETERS

PHYSICAL ANALYSIS		Name of the Method and the Standard	TS 266 April-05
pH		SM (2005) 4500 - H+ .B. Electrometric Method	6,5-9,5
Turbidity	NTU	SM (2005) Turbidity 2130 B. Nephelometric Method	1
CHEMICAL ANALYSIS			
Ammonia	mg/L	S.M 4500-NH3-D Ammonia-Selective Electrode Method	0,5
Ammonium	mg/L	TS EN ISO 14911 Nisan 2000 - Ion Chromatography	0,5
Bromate	mg/L	SM (2005) 4110 B-Ion Chromatography with Chemical Suppression of Eluant Conductivity	0,01
Bromide	mg/L	SM (2005) 4110 B-Ion Chromatography with Chemical Suppression of Eluant Conductivity	
Fluoride	mg/L	SM (2005) 4110 B-Ion Chromatography with Chemical Suppression of Eluant Conductivity	1,5
Phosphate	mg/L	SM (2005) 4110 B-Ion Chromatography with Chemical Suppression of Eluant Conductivity	
Calcium	mg/L	TS EN ISO 14911 Nisan 2000-Ion Chromatography	
Chloride	mg/L	SM (2005) 4110 B-Ion Chromatography with Chemical Suppression of Eluant Conductivity	250
Lithium	mg/L	TS EN ISO 14911 Nisan 2000-Ion Chromatography	
Magnesium	mg/L	TS EN ISO 14911 Nisan 2000 - Ion Chromatography	
Nitrate	mg/L	SM (2005) 4110 B-Ion Chromatography with Chemical Suppression of Eluant Conductivity	50
Nitrite	mg/L	SM (2005) 4110 B-Ion Chromatography with Chemical Suppression of Eluant Conductivity	0,5
Potassium	mg/L	TS EN ISO 14911 Nisan 2000-Ion Chromatography	
Sodium	mg/L	TS EN ISO 14911 Nisan 2000-Ion Chromatography	200
Sulfate	mg/L	SM (2005) 4110 B-Ion Chromatography with Chemical Suppression of Eluant Conductivity	250
Free Chlorine	mg/L	SM (2005) 4500-Cl F DPD Ferrous Titrimetric Method	
BACTERIOLOGICAL ANALYSIS			
T.Coli.Bact.	in 100mL	SM (2005) 9222 B Total Coliform Membrane Filter Procedure	0
E coli	in 100mL	SM (2005) 9222 B Total Coliform Membrane Filter Procedure	0
T.Coli.Bact.	EMS/100mL	SM(2005) 9221 B-C-E Total Coliform Fermentation Technique	0
Fecal Coliform	EMS/100mL	SM(2005) 9221 B-C-E Total Coliform Fermentation Technique	0
Total Bacteria Count (HPC) at 35 °C	kob/mL	SM(2005) 9215 A-B Pour Plate Method	-
Total Bacteria Count (HPC) at 22 °C	kob/mL	SM(2005) 9215 A-B Pour Plate Method	-

CHAPTER 5

WATER QUALITY CONTROL LABORATORY HAS BEEN ACCREDITED

Our laboratory has achieved certificate of compliance to "TS EN ISO/IEC 17025:2005 General Conditions for the Competency of Experiment and Calibration Laboratories" in 2007 and proved its international service quality.

**PARAMETERS AND METHODS USED IN THE ANALYSES IN WATER QUALITY CONTROL LABORATORY,
GENERAL DIRECTORATE OF İSKİ**

Order No	Name of the analysis	Analysis method Used	Mane of the Used Standard
1	*pH	pH Measurement Method	SM [2005] H+ 4500 - H.B. Electrometric Method
2	*Turbidity	Nephelometric Method	SM [2005] Turbidity 2130 B. Nephelometric Method
3	* Free chlorine	DPD and FAS Titrimetric Method	SM [2005] 4500-Cl F DPD Ferrous Titrimetric Method
4	Free chlorine	Visual Comparison Technique	EPA- Water-quality-sampling-manual Part D 8.5
5	Total chlorine	Visual Comparison Technique	EPA- Water-quality-sampling-manual Part D 8.5
6	Color	Color Analysis Method	SM [2005] 2120 B. Visual Comparison Method
7	*Ammonia	Ion Selective Method	SM [2005] 4500- NH3 D Ammonia-Selective Electrode Method
8	Nitrate	Ion Selective Method	SM [2005] 4500 NO3- D Nitrate Electrode Method
9	Total dissolved Solids	Electrometric Method	SM [2005] 1030-E Checking Correctness of Analysis
10	Alkalinity	Titrimetric Method	SM [2005] Alkalinity 2320 B. Titrimetric Method
11	*Ammonium	Cation Analysis Method with Ion Chromatography Equipment	TS EN ISO 14911 April 2000
12	*Bromate	Anion Analysis Method with Ion Chromatography Equipment	SM [2005] 4110 B-Ion Chromatography with Chemical Suppression of Eluant Conductivity
13	*Bromide	Anion Analysis Method with Ion Chromatography Equipment	SM [2005] 4110 B-Ion Chromatography with Chemical Suppression of Eluant Conductivity
14	*Fluoride	Anion Analysis Method with Ion Chromatography Equipment	SM [2005] 4110 B-Ion Chromatography with Chemical Suppression of Eluant Conductivity
15	*Phosphate	Anion Analysis Method with Ion Chromatography Equipment	SM [2005] 4110 B-Ion Chromatography with Chemical Suppression of Eluant Conductivity
16	*Calcium	Cation Analysis Method with Ion Chromatography Equipment	TS EN ISO 14911 April 2000
17	*Chloride	Anion Analysis Method with Ion Chromatography Equipment	SM [2005] 4110 B-Ion Chromatography with Chemical Suppression of Eluant Conductivity
18	*Lithium	Cation Analysis Method with Ion Chromatography Equipment	TS EN ISO 14911 April 2000
19	*Magnesium	Anion Analysis Method with Ion Chromatography Equipment	TS EN ISO 14911 April 2000
20	*Nitrate	Anion Analysis Method with Ion Chromatography	SM [2005] 4110 B-Ion Chromatography with hemical Suppression of Eluant Conductivity
21	*Nitrite	Anion Analysis Method with Ion Chromatography Equipment	SM [2005] 4110 B-Ion Chromatography with hemical Suppression of Eluant Conductivity
22	*Potassium	Cation Analysis Method with Ion Chromatography Equipment	TS EN ISO 14911 April 2000
23	*Sodium	Cation Analysis Method with Ion Chromatography Equipment	TS EN ISO 14911 April 2000
24	*Sulphate	Anion Analysis Method with Ion Chromatography Equipment	SM [2005] 4110 B-Ion Chromatography with Chemical Suppression of Eluant Conductivity
25	Conductivity	Conductivity method	SM [2005] Conductivity 2510 B. Laboratory Method
26	*Total hardness	Titrimetric Method	SM [2005] 2340 C. EDTA Titrimetric Method
27	Geosmin	Mikroextraction Method	SM [2005] 6040 D- Solid Phase Microextraction (SPME)

28	MI B	Mikroextraction Method	SM [2005] 6040 D- Solid Phase Microextraction [SPME]
29	*Total organic carbon	High- Temperature Combustion Method	SM [2005] 5310 B. High- Temperature Combustion Method
30	Permanganate index	Titrimetric Method	TS 6288 EN ISO 8467/Nisan 1998
31	Sulphate	Turbidimetric Method	SM [2005] 4500 S04-2 E. Turbidimetric Method
32	Calcium	Titrimetric Method	SM [2005] 3500 Ca B. EDTA Titrimetric Method
33	Magnesium	Calculation Method	SM [2005] 3500 Mg B. Calculation Method
34	Chloride	Titrimetric Method	SM [2005] 4500-Cl- B. Argentometric Method
35	*UV254-Absorbing Organic Constituents	Spectrophotometric Method	SM [2005] 5910-B Ultraviolet Absorption Method
36	Cyanide	Ion Selective Electrode Method	SM [2005] 4500 CN- -F. Cyanide -Selective Electrode Method
37	Fluoride	Ion Selective Electrode	SM [2005] 4500 4500 F- C. Ion- Selective Electrode Method
38	Chloroform	Purge and Trap Gas Chromatographic Method	SM [2005] 6232 D. Purge and Trap Gas Chromatographic Method
39	Dichlorobromomethane	Purge and Trap Gas Chromatographic Method	SM [2005] 6232 D. Purge and Trap Gas Chromatographic Method
40	Dibromochloromethane	Purge and Trap Gas Chromatographic Method	SM [2005] 6232 D. Purge and Trap Gas Chromatographic Method
41	Bromoform	Purge and Trap Gas Chromatographic Method	SM [2005] 6232 D. Purge and Trap Gas Chromatographic Method
42	Determination of Trihalomethanes (THM's)	Purge and Trap Gas Chromatographic Method	SM [2005] 6232 D. Purge and Trap Gas Chromatographic Method
43	Algae toxins (LR-RR)	Determination by HPLC with UV-VIS (Ultraviolet-Visible) detectors	Analyst, July 1994, V.119 (1525-1529) HPLC App.No:H10
44	*Total coliform	Membrane Filter Method	SM [2005] 9222 B Standard Total Coliform Membrane Filter Procedure - TS EN ISO 9308-1:2000
45	*E.Coli	Membrane Filter Method	SM [2005] 9222 B Standard Total Coliform Membrane Filter Procedure - TS EN ISO 9308-1:2000
46	*Total coliform	Multi-Tube Method	SM[2005] 9221 B-C- Standard Total Coliform Fermentation Technique-SM[2005] 9221 E. Fecal Coliform Procedure
47	* Fecal Coliform	Multi-Tube Method	SM[2005] 9221 B-C- Standard Total Coliform Fermentation Technique-SM[2005] 9221 E. Fecal Coliform Procedure
48	*Total bacteria	Pour Plate Method	SM[2005] 9215 A-B Pour Plate Method
49	Chlorophyll (a)	Acetone Extraction Method with Spectrophotometer	SM [2005] 10200 H2 Chlorophyll
50	Giardia	Determination Method under Fluorescence Microscope with IMS (Immunomagnetic Separation) and FA	SM[2005] 9711 B Giardia and Cryptosporidium Methods EPA-1622-1623
51	Cryptosporidium	Determination Method under Fluorescence Microscope with IMS (Immunomagnetic Separation) and FA	SM[2005] 9711 B Giardia and Cryptosporidium Methods EPA-1622-1623
52	Phytoplankton Count and Species Identification Method	Phytoplankton Counting Method with Sedgwick-Rafter(S-R) counting cell	SM [2005] 10200 F- Phytoplankton Counting Techniques
53	*Determination and Counting of Intestinal Enterococci by Membrane Filter Method	Membrane Filter Method	TS EN ISO 7899-2-Detection and enumeration of intestinal Enterococci - Membrane Filter Method
54	Salmonella	Membrane Filter Method	TS ISO 6340/ April 1999 Search for Salmonella

2007 WATER QUALITY REPORT OF İSTANBUL 2007

Parameter	TURKISH STANDARTS TS 266 April 05	WORLD HEALTH ORGANIZATION (WHO) 1999	USA ENVIRONMENTAL PROTECTION AGENCY (EPA) 2002	EUROPEAN COMMUNITY (EC) 1998	MEAN QUALITY VALUES OF REFINING FACILITY (GENERAL OUTPUT)						
					Büyükçekmece	İkitelli	Kağıthane	Ömerli (Emirli)	Ömerli (Mur.-Orhaniye)	Ömerli (Osmaniye)	Elmalı
Turbidity	1	5	1	1	0,2	0,3	0,1	0,3	0,4	0,5	0,3
PRIMARY STANDARDS (MICROBIOLOGICAL), EMS/100 ml											
Coliform Bacteria	0	0	0	0	0	0	0	0	0	0	0
PRIMARY STANDARDS (Disinfection by-products), mg/l											
Total Trihalomethane	-	460	100	100	45	43	35	24	14	9	27
Bromate	10	25	10	10	10	7	14	7	5	5	4
PRIMARY STANDARDS (INORGANIC CHEMICALS), mg/l											
Aluminum	0,20	0,20	0,20	0,20	0,043	0,044	0,034	0,000	0,066	0,074	0,070
Arsenic	0,01	0,05	0,01	0,01	0,000	0,000	0,000	0,000	0,000	0,000	0,000
Barium	-	0,7	2	-	0,048	0,028	0,028	0,039	0,025	0,027	0,027
Cadmium	0,005	0,005	0,005	0,005	0,000	0,000	0,000	0,000	0,000	0,000	0,000
Chromium (Total)	0,05	0,05	1	0,05	0,000	0,000	0,000	0,000	0,000	0,000	0,000
Fluoride	1,5	1,5	2	1,5	0,21	0,10	0,10	0,10	0,07	0,07	0,12
Lead	0,01	0,05	0,015	0,01	0,000	0,000	0,000	0,000	0,000	0,000	0,000
Mercury	0,001	0,001	0,002	0,001	0,000	0,000	0,000	0,000	0,000	0,000	0,000
Nitrate (NO3-)	50	50	45	50	3,5	2,8	2,6	4	4,5	4	5
Selenium	0,01	0,01	0,05	0,01	0,000	0,000	0,000	0,000	0,000	0,000	0,000
Silver	-	--	0,1	-	0,000	0,000	0,000	0,000	0,000	0,000	0,000
Antimony	0,005	0,005	0,006	0,005	0,000	0,000	0,000	0,000	0,000	0,000	0,000
Beryllium	--	--	0,004	-	0,000	0,000	0,000	0,000	0,000	0,000	0,000
PRIMARY STANDARDS (RADIOLOGICAL), pCi/l											
Gross Alfa	0,1	-	-	-	-	-	-	-	-	-	-
Gross Beta	1	-	-	-	-	-	-	-	-	-	-
SECONDARY STANDARDS (AESTHETIC), mg/l											
Chloride	250	250	250	250	91	113	168	40	25	24	46
Color (unit)	20	15	15	-	2,6	2,7	2,5	2,5	2,7	2,6	2,6
Copper	2	--	1	2	0,001	0,001	0,001	0,000	0,001	0,002	0,001
Detergents	-	--	0,5	-	-	-	-	-	-	-	-
Iron	0,2	--	0,3	0,2	0,003	0,003	0,004	0,005	0,005	0,005	0,007
Manganese	0,05	0,5	0,05	0,05	0,000	0,000	0,000	0,000	0,01	0,01	0,01
Odor producing Geosmin	-	-	-	-	-	-	-	-	-	-	-
Substances MIB	-	-	-	-	-	-	-	-	-	-	-
PH	6,5-9,5	6,5 - 8,5	6,5 - 8,5	6,5 - 9,5	7,2	7,2	7,1	7	6,8	6,9	7
Sulphate	250	250	250	-	83	50	58	13	32	33	58
Total Soluble Substance	-	1000	500	-	304	325	392	182	135	141	204
Zinc	-	-	5	-	0,001	0,004	0,007	0,037	0,005	0,003	0,007
ADDITIONAL PARAMETERS, mg/l											
Calcium	-	--	-	-	48	48	48	53	33	36	42
Hardness (as CaCO3)	-	500	-	-	177	158	177	167	110	115	141
Magnesium	-	-	-	-	13,8	10,3	12,7	6	5,3	5,5	9,3
Potassium	-	-	-	-	4,8	4,3	4,8	1,5	1,8	1,8	3
Sodium	200	200	-	200	51,2	64,9	86,8	12,9	12,9	12,9	27,1
Free Chlorine (3)	-	5	4	-	1,3	1,3	1,1	1,1	1,2	1,1	1,2
Ammonium	0,5	1,5	-	0,5	< 0,03	< 0,03	< 0,03	< 0,03	< 0,03	< 0,03	< 0,03

1- Turbidity, pH, Chloride, Total Hardness and Residual Chloride values are the monthly mean values of the Operation Laboratory of the Facility. 2- Other Parameters are the monthly mean values of the analyses performed in Waste Water and Water Quality Laboratories. 3- Chloride value shown here is the facility output value and free chloride decreases with time. This value is chosen to have free chloride in the utmost point of the system, Free chloride values are generally measured in within the range of 0,1-0,5.

Anion-Cation analysis

Chromatography is the general name of the methods that enable the separation of the components of a mixture with the help of a phase and therefore yielding qualitative and quantitative analyses. In ICS 3000 ion chromatography device, anions (F, Cl, NO₂, Br, BrO₃, NO₃, SO₄, PO₄) and cations (Li, Na, NH₄, K, Mg, Ca) can be analyzed simultaneously.

This method has some advantages over wet chemical analyses. These are;

- Both anions and cations can be analyzed within half hour with a single injection.
- The sample is injected directly to the column.
- Can analyze samples over wide study intervals
- Obtaining sensitivity as the level of $\mu\text{g/l}$
- Performing analyses even with low sample volumes such as 15-150 μL .

Ion chromatography device

Trihalomethanes (THM's)

Analysis of trihalomethanes (THM's), which are among the primary parameters of drinking water as a by-product of disinfection of water is performed on Gas chromatography.

The purpose of this analysis performed on this device is following the Trihalomethanes (Chloroform, Dichlorobromomethane, Dibromochloromethane, bromoform) concentration in water and determine the possible increases and informing the drinking water purification facility and having them take the necessary precautions.

Total Organic Carbon (TOC)

Total organic carbon analysis is analysis that determines the complete organic substances in water samples. Total organic carbon quantities' being high makes disinfection difficult, and at the same time increases the formation of trihalomethane.

Bromate

In bromide containing waters, it is by-product that is formed as a result of ozoning of waters for disinfection purposes. It has cancerogenic effect and should not be present in drinking waters above the level of $25\mu\text{g/L}$.

Lithium

One of the minerals that water contains as a trace element ($0,2\text{ mg/L}$ or less).

Taste and Odour Parameters (MIB-Geosmin)

In superficial water sources, with the increasing temperature, algae concentration also increases. Some algae species may give the water a foul smell and odor.

The purpose of the analysis performed on this device is following the concentration of 2-Methylisoborneol and Geosmin that lead to foul smell and odor in water and determining the possible increases and informing the drinking water purification facility and having them take the necessary precautions (such as the addition of powdered active carbon).

Gas Chromatography - Mass spectrometry (GC-MS) system

Gas Chromatography (GC) system

**Nitrite**

Drinking water and potable water normally should not contain nitrite. Therefore, nitrite analysis in waters may be measured qualitatively (present-absent). Ammonia, the first decomposition product of organic substances is oxidized and then nitrites form. They turn into nitrates in a very short time. Though hazardous for human health, effects of nitrite in trace amounts are not yet known.

Sodium

Excessive amounts lead to taste problem. Leads to the feeling of salinity. Excess sodium may originate from household and industrial pollution, soil structure and sea additives.

Potassium

As in sodium, results in taste problems. Originates from industrial pollution, agricultural manures and soil structure.

Phosphate

Has the potential to form stones. Presence of phosphorus in lakes is vital for the life in lakes, however, excessive phosphorus leads to the over reproduction of algae in lakes. Presence of phosphorus draws attention to other parameters as well. It is a parameter that can be monitored. Originates from household wastes, manures and industrial wastes.

Calcium

Calcium does not have a direct harmful effect on the body. It is even claimed to be beneficial for bone structure. However, it constitutes a problem regarding drinking. On the other hand, stone formation potential of the water increases. Because, if it is found in low quantities, this may result in abrasive effect. Calcium in water is due to the structure of the soil that water passes through.

Magnesium

When present in drinking water in excessive amounts, results in damages in the eyes. Diarrheic effect is also seen. Magnesium in

water is due to the structure of the soil that water passes through. Gives bitterness to the water.

Fluoride

Fluoride present in waters may be beneficial or hazardous according to its quantity. The recommended amount for drinking water is 1 mg/L. This value is known to be beneficial for teeth and decrease tooth decay. In addition, the damages that high levels of fluoride containing waters cause in humans are proven with researches. For example, in a study conducted in children under the age of 9 demonstrated that water containing 2 mg/L fluoride leaves brown stains on teeth and water containing 4 mg/L fluoride leads to bone disorders. In this case, according to the results of this research, waters containing more than 1 mg/L fluoride should be purified. fluoride should be purified.

Total Hardness

Hardness is the result of the presence of (+2) ions (Ca^{+2} , Mg^{+2} , Sr^{+2} , Fe^{+2} , Mn^{+2} etc.) that dissolve in water. Since Ca^{+2} and Mg^{+2} ions are found in natural waters more than the other ions, hardness is mainly expressed as the Sum of Ca^{+2} and Mg^{+2} ion concentrations. Since other ions are generally in complex form, they do not contribute that much to hardness.

Hardness of water is the measurement of its soap precipitation capacity. Soap is commonly precipitated by calcium and magnesium ions in water. Sabun suda yaygın olarak kalsiyum ve magnezyum iyonları ile çöktirilir. Though ions of some other metals (Al , Fe , Mn , Sr , Zn) precipitate soap as well, since they are usually in the form of a complex, they do not contribute that much to hardness.

Waters may be classified according to their hardness degrees as follows.

Total Hardness (mg CaCO_3/l)	Rating
0-75	Mild water
75-100	Water of moderate hardness
100-300	Hard water
$\rightarrow 300$	Very hard water

consideration that organic substances may also be present in the water, it should be considered as an element of contamination. It is difficult to detect ammonia in chlorinated waters. Because, ammonia may be present as chloramine compounds. Therefore, before analyzing chlorinated waters, dechlorination should be performed. If ammonia is present in water, this shows that water is contaminated by organic substances at a close point. Ammonia is formed in the first phase of the cycle where organic substances pass through nitrification and become completely harmless. If the period after the organic substance is mixed with water is short, i.e., not enough to complete nitrification, ammonia is detected in the water.

Ammonium

Constitutes taste and odor problem; has a negative effect on human health. Presence of ammonium in water shows that water is contaminated by household wastes. Therefore, these kinds of waters are potentially diseases forming. Originates from household and industrial pollution, manures.

Nitrate

These are the last production of oxidization of nitrous organic compounds. Nitrate is usually found more in well waters. Especially leads to a disease called blue-baby in babies. This disease may even lead to death in babies whose bodies start to turn into purple. Nitrates may come from soil, however if it originates from ammonia and nitrite, precautions should be taken. Because the presence of nitrites means that the water is contaminated. If nitrites are present with large amounts of organic substances, then a bigger contamination

is in question. Ammonia leads to reproduction of some bacteria species and these bacteria give the water foul smell.

Nitrogen types are given to the recipient environment in excess amounts, they are used by the organisms. This leads to eutrophication (lack of oxygen due to algae boost) in this recipient environment.

According to (TS-266) 1997 standards, nitrate is recommended to be present in drinking water below 25 mg/L. This value should not exceed 50 mg/L.

Chloride

Chloride is the type of ion that is very commonly present in all natural or used waters. May be added to water as a result of dissolving from underground formations or due to interventions of salt water-fresh water. 6 g chloride/day on the average is excreted via human urine per person.

does not lead to a health problem at normal concentrations.

However, at concentrations above 250 mg/L, salty taste forms. Since chloride increases the conductivity of water, it facilitates corrosion. High chloride concentration manifests its presence in water with the adverse effect on taste, corrosive tendency or softening process.

Ammonia

Normally, ammonia should not be present in drinking water. If ammonia is present in water, this shows that water is contaminated by organic substances at a close point. Ammonia is formed in the first phase of the cycle where organic substances pass through nitrification and become completely harmless. If the period after the organic substance is mixed with water is short, i.e., not enough to complete nitrification, ammonia is detected in the water. In addition, ammonia may be found in the water due to the denitrification of nitrites and nitrates in water with bacteria. Without examining the presence of other chemical substances, it can not be decided that the water detected to contain ammonia is contaminated. Ammonia does not have hazardous effects on health alone. However, since it leads to the

Chemical Analyses

Free Chloride

Chlorination of drinking water is a disinfection process performed to eliminate pathogenic microorganisms that water contains and make water safe for drinking. According to the character of raw water, dosing is performed in the range of 0,5 ppm and 5 ppm.

Sulphate

Sulphate is one of the most important ions that mix with surrounding waters naturally. All natural waters contain sulphate in varying amounts. Sulphate content of some industrial waste waters is high and when they mix with natural waters, they increase these waters' sulphate contents as well. Sulphur (S (-II)) compounds are important contaminants due to problems they establish like taste, odor, toxicity and corrosion with various reactions. The meaning of high sulphate content of water means high hardness, high sodium salt and high acidity.

Since sodium sulphate and magnesium sulphate show diarrheic effect in humans, it is restricted to an upper limit of 250 mg/L. For animals, this limit is determined to be 1000 mg/L. In addition, sulphates give a slight bitter taste to the water.

Conductivity

Conductivity is a function of ions dissolved in water. Therefore, it is a monitoring parameter. Electrical conductivity of water. Waters that do not contain inorganic substances (for example, deionized pure water) have difficulty in conducting electricity. Conductivity also decreases with increasing temperature. Therefore, measurement of electrical permeability of water may be considered as a measurement of the electrolyte amount in water. The measurement is performed between 18-20°C.

Conductivity Device

TDS (total dissolved solids)

Total dissolved solids in water contain inorganic salts and organic substances to a lesser degree. Both surface and underground waters dissolve minerals from soils and stone materials that they come in contact with. Dissolved inorganic materials are present in waters as ions. The most recognized general ions in water are as follows;
Cations: Ca^{+2} , Mg^{+2} , Na^{+2} , Fe^{+2} , Mn^{+2}
Anions: HCO_3^- , Cl^- , SO_4^{+2} , NO_3^- , CO_3^{+2}

In addition, waters may contain heavy metal ions (lead, mercury, cadmium etc.) and organic substances. Dissolved organic chemicals like (pesticides, herbicides) have toxic effect on humans and animals even in low quantities. Most of the dissolved organic substances like trihalomethanes (THM) and dioxine lead to cancer. These kind of organic substances are in the form of ions and have a low concentration. Above mentioned ions provide electrical conductivity in water. It is corrosive for metallic surfaces in high quantities. Excess TDS may lead to stratification in pipes; high concentration in drinking water may have diarrheic effect.

Total dissolved solids may be removed only by reverse osmosis and demineralization processes.

Turbidimeter

Turbidity

Turbidity is the measure of light permeability of waters containing suspended solid materials. Turbidity may be composed of many things, from suspended substances in the water to apparent large precipitates. Substances like sand, clay, silica, calcium carbonate, iron, manganese, sulphur etc. causes turbidity. Turbidity, especially abundant in inland waters is based on soil transported by rains and household-industrial waste waters that mix with inland waters. In addition, during this contamination, inorganic substances are also mixed into the water as organic substances. The presence of these substances supports the bacteria growth in water. And bacteria formation increases the turbidity of water. For example, algae, using substances like N.P grow and cause turbidity in water. In addition, increasing temperature also increases the microorganism activities. As a consequence, the reason of turbidity may completely be based on inorganic substances or be based on and many organic substances present in the nature.

Color

Color increase in water is not desired. Because this situation shows the presence of metals like iron, manganese, chromium, nickel and organic compounds that have dissolved in water. Color is also undesired in respect to aesthetics. Colored water may result from decomposition of colored plants in water, the structure of the soil, household and industrial pollution.

Color Device

STUDIES OF THE CHEMISTRY LABORATORY

Physical Analyses

Ph

PH is a logarithmic measure of acidity and alkalinity of water. Expresses the H^+ ion concentration present in the solution.

if $pH < 7$, then the medium is acidic.

if $pH > 7$, then the medium is basic.

pH measurement is one of the most frequently used and most important parameters that is used in the analysis of water chemistry. It is based on CO_3^{2-} , HCO_3^- , CO_2 , and OH^- ions. In waters having low pH, i.e. acidic, CO_2 and HCO_3^- ions are dominant. Water having a pH lower than 6.5 is acidic and thus has corrosive property. Due to this property, has corrosive effect on the metals present in purification facilities, network system and houses. In water, havin a pH higher than 9.5 itaste and odor problems arise; it gives the water a soapy, slippery sensation. In addition, these waters have the property of producing stones. Drinking and potable water are required to have

pH meter

HPLC Device

Algae Toxins Analysis Method

Some species of Blue-Green algae produces toxins. Increase in toxin concentration during algae boost (overproduction) causes diseases in humans and animals. In our laboratory, Microcystin LR and RR analysis, which are among the most important algae toxins is performed.

Method: Describes the algae toxin determination method in the water sample by HPLC.

Based on determination with DAD (Diode Array Dedector) and UV-VIS (Ultraviolet-Visible) detector by passing Algae toxins through HPLC column.

Source: Analyst, July 1994, V.119 [1525-1529] HPLC App. No:H10
Blue-Green Algae are the most widespread algae species in drinking

Phytoplankton Counting and Species Determination Method

water. High levels of algae containing water obstructs the fFilters and makes the facility work with reduced capacity. Causes foul smell, taste and visual pollution.

Method: Describes the method of microscopical Phytoplankton count and species identification

Source: SM (2005) Sec. 10-13 10200 F- Phytoplankton Counting Techniques

Fluorescence Microscope

Method: Determination and Counting of Intestinal Enterococci with Membrane Filter Method gives the direct count of microorganisms in the sample.

Source: TS EN ISO 7899-2/April 2002

Giardia and Cryptosporidium Study Method

Among the Pathogenic protozoa; Giardia cysts and Cryptosporidium oocysts are resistant to chemical disinfectants compared to bacteria and viruses. Cysts are distributed by the feces of infectious persons. Cysts are transmitted by oral ingestion of contaminated water and food. Causes diarrhea, gastroenteritis and absorption disorders. Method: Determination of Giardia cyst and Cryptosporidium oocysts in Samples with Filtration, IMS (Immunomagnetic separation) and FA (Fluorescence antibody) methods:

Source: -SM(2005) Sec.9-162 9711 B Giardia and Cryptosporidium Methods

IMS System

Computer assisted
Fluorescence microscope

Total Bacteria Count Analysis method

With Total Bacteria, only microorganism group that can develop under the specified conditions are meant.

It covers the Total count of all colonies, reproducing on the surface of agar as a result of its incubation at 35(0.5°C and 22)0.5°C. When total bacteria count increases, necessary precautions should be taken.

Method: Pour Plate Method - Accredited-Is a method that gives the number of Total Mesophilic Aerobic Bacteria in the sample

Source: SM (2005) Sec. 9-37 9215 A-B Pour Plate Method

Determination of Intestinal Enterococci by Membrane Filter Method

Enterococci are present in the intestines of warm blooded animals including humans; they are more resistant to environmental conditions than E.coli and Coliform bacteria. Since they are more resistant, they are important in the determination of the efficiency of the purification process. In respect to water analysis, enterococci are considered to be the indicator of fecal contamination. Their presence in Drinking/Potable water is the indicator of direct or indirect contamination of the sample by feces. Usage of these waters, contaminated by bacteria which are present in normal intestinal flora results in gastrointestinal diseases.

Distribution of Mediums by Dosimeter

Autoclave

Performance of Multi-Tube Culture-Passage Performance

Total Coliform and Fecal Coliform Detection with Multi Tube Method

Fecal Coliform, is a subgroup of Coliform group and originates from feces. Bacteria defined as Fecal coliforms are firstly E.coli. Its presence in drinking/potable water is the indicator of direct or indirect contamination of the sample with feces.

Usage of water contaminated by these bacteria, which are present in normal intestinal flora results in gastrointestinal diseases.

Method: Multi tube Fermentation Technique - Accredited - Bacteria count is determined by using the Tables as microorganism intensity in sample MPN (Most Probable Number) and Turkish EMS (En Muhtemel Sayı).

Source: SM (2005) Sec. 9-49 9221 B-C- Standard Total Coliform Fermentation Technique

-SM (2005) Sec 9-56 9221 E. Fecal Coliform Procedure

STUDIES OF THE BACTERIOLOGY LABORATORY

Total Coliform and E. Coli Determination by Membran FFilter Method

Since the determination and detection of pathogenic bacteria in water is difficult, instead, coliform bacteria are used as indicator microorganisms since they can be easily identified and counted in water according to International Standards (ISO, EPA, WHO). E.coli is considered to be proof of contamination of water with feces. Since these bacteria are present in water, it is also considered that pathogenic bacteria having a high risk of disease would also be present. Usage of water contaminated by these bacteria, which are present in normal intestinal flora results in gastrointestinal diseases.

Method: Membrane FFilter Method-Accredited- It is a method that directly determines the organisms present in the sample.
Source: -SM (2005) Sec. 9-60 9222 B Standart Total Coliform Membrane FFilter Procedure
- TS EN ISO 9308-1:2000

Membrane Fiter filtration
system in operation

MiniAPI Identification
system used for the
Identification of bacteria

Important Precautions in Sample Collection

General collection rules that you have to abide by while collecting the samples that you will send to our laboratory for analysis are as follows. In order to obtain accurate results, these rules have to be strictly followed.

Sample Labelling: Sample type has to be specified on the sample containers. Descriptive information about the sample (Sample name, institution name, sample code, sampling date, name of the collector) has to be written on the sample container.

Sample Quantities and Sample Collection Containers: Water samples have to be collected in clean plastic barrels, glass jars having a volume of at least 1 liter. If possible, the cap should be closed so as not to contain air between the sample and the cap.

Sample Containers

Bacteriological sample collection;

- Tap should be turned on and water should flow for 1-2 minutes
 - With a cotton, wrapped around a metal stick and dampened with alcohol, or with a lighter, outlet of the tap is burnt (Picture 1)
 - Tap is turned on again
 - In order to prevent contamination, to avoid touching the cap and mouth, the container is held at the middle part, and is filled to the bottle neck with water to provide good mixing. (Picture-2/3)
 - Water sample should be delivered to the laboratory within 6 hours and analysis should be started.
 - Under extraordinary conditions, water sample can be delivered to the laboratory in a fridge within 24 hours.
 - Disposable, thiosulphate ($\text{Na}_2\text{S}_2\text{O}_3$) containing propylene sample containers are used for the collection of samples for bacteriological analysis (thiosulphate is used to prevent the negative effect of chloride
- Source: TS 5106-ISO 5667-3/April 1997 on bacteria

Picture-1

Picture-2

Picture-3

Source: TS 5106-
ISO 5667-3/Nisan 1997

CHAPTER 4

LABORATORY STUDIES

One of the fundamental elements of a healthy life is having clean drinking and potable water. It can be decided by BACTERIOLOGICAL AND CHEMICAL analyses whether water is healthy, clean and potable.

Sample Admission Unit

ISKİ personnel collect water samples from treatment plants, lakes, dams, springs and wells providing water, and also from the city's water network within the borders of Municipality and from water depot outlets and network outlets of schools and hospitals on a daily, weekly, monthly basis according to "TSE 5106 Standards of Sample Collection, Sample Transportation and Storage Conditions", with a pre-determined schedule, and deliver all samples to the laboratory

under appropriate conditions. The delivery of the samples brought to the laboratory is taken by the Sample Admission and Extermination Officer, by controlling their compatibility to ISKİ Water and Waste Water Sample Quantities Storage Conditions and Periods.

These technical data which are related to the samples having an important form of delivery from the customer is presented in the website for informative purposes.

Sample Admission
Room

CHAPTER 3

PROPERTIES OF DRINKING WATER

Water has to be odourless, colorless and clear and provide a pleasant drinking experience.

Materials like phenols and oils which give water a foul taste and odour should not be present in water. Water has to be clear and colorless. Micro-organisms causing diseases should not be present in water. Micro-organisms like vibrio cholera, salmonella typhi, hepatitis virus in water causes diseases. Drinking waters should absolutely be free of contamination.

Hazardous chemicals should not be present in water. Some chemical substances might have poisonous affect. Like arsenic, cadmium, chromium, lead, mercury. In addition, substances like barium, nitrate, fluoride, radioactive substances, ammonium, and chloride have negative effects on health if they are present above limit values. At the same time, some are the indicators of contamination of water with dirty waters.

Water has to be suitable for its purpose of use.

Iron, manganese presence and hardness of water are important in drinking water and in potable water used in the industry. Water should not be aggressive.

Aggressiveness of water is caused by the equilibrium of free carbon dioxide (CO_2) and bicarbonate (HCO_3^-) ion. Aggressiveness of water leads to the corrosion of water pipes. In addition, if pipes are corroded, then elements released from the pipe corrupt the quality of water.

daily, weekly and monthly periods, performing chemical and bacteriological analyses of these samples; reportating, inspecting whether they are in accordance with TS-266 "Drinking and Potable Water Standards" and performing further necessary analyses. In addition, water depot outlets and network outlets of schools and hospitals are controlled and District Educational Directorates and District Public Health Directorates are alerted accordingly by İSKİ and disinfection is also performed, as a result, İSKİ works in coordination with other local institutions in order to provide standardized water to the network.

Tap water of İstanbul is analyzed in modern laboratories, by obtaining samples from 350-400 different points on a daily basis. Analysis results demonstrate that tap water of İstanbul is above criteria of Turkish Standardization Institution (TSE 266), European Community (EC), World Health Organization (WHO) and U.S. Environmental Protection Agency (EPA). Today, with these studies of İSKİ, healthy water that can safely be drank flows from every tap of every house in İstanbul.

Picture-1

cleaned. After all these processes, water is chlorinated and distributed by clean water network system.

The most important purpose in water purification is the protection of human health. With this purpose, İSKİ uses state of the art technology to purify water at its drinking water treatment plants having a universal large production capacity.

Quality changes in raw water and with the help of parameters obtained in every step of purification, output values of clean water (i.e. pH, turbidity, ammonia, color, chlorine, alkalinity, conductivity, hardness etc.) are continuously monitored in order to obtain standardized clean water, and necessary changes in the purification facility are made. All factors that can affect the quality of water are handled and controlled starting from the source itself.

Picture-1

İSKİ analyses water in its laboratories equipped with the state of the art devices with scientifically valid analysis methods and reports according to national and international standards.

İSKİ is responsible for collecting water samples from the Treatment Plants, Lakes, Dams, Springs and Wells providing water, and also from city water network within the borders of the Municipality with

Kağıthane
Water treatment facility

CHAPTER 2

THE ADVENTURE OF WATER POURING FROM FOUNTAINS

Water, that has an important place in our lives, is one of the indispensable elements of a healthy life. Natural spring waters are one the basic requirements of the human body with the minerals they contain. They play a very important role in the nourishment of humans; energy production; transportation, absorption and metabolism of nutrients; waste disposal; transportation of oxygen to the cells and removal of carbon dioxide. Water does not only consist of water, it also contains minerals like calcium, magnesium, and sodium. The consumption of healthy water having a high level of minerals is important for a healthy life of good quality in all ages. Treatment of drinking water

Water taken from the spring and brought to the purification facility is firstly passed through aeration pools. Here, dissolved oxygen ratio in water is increased, undesirable substances leading to foul taste and odor are removed and iron oxidation is performed. Then, contaminants are oxidized in the ozonization unit. The next step is fast and slow mixing pools. Here, suspended solid materials that can not be precipitated with the addition of chemical substance are turned into flocks and eliminated in precipitation pools. Organic and inorganic contaminants that could not be removed are passed through fast sand filters, the next unit in the process, thus water is clarified and

CHAPTER 1

THE ADDRESS OF WATER AND QUALITY: İSKİ

Great civilizations have always been established near waterfront, and dominated by water. A city, existing via water has been the centre of civilizations. Water has brought the city of İstanbul such a civilization that it has made it the capital of water civilization.

With the increase in water consumption due to urbanization, all advanced civilizations made great efforts to bring water to the cities, to store, refine and distribute water that has been the most important element in the life of the city and people.

Romans established an advanced and successful water system by bringing water they have obtained from sources to various fountains via aqueducts. During Byzantine period, usage of cisterns in which water was stored was more common. Ottoman sovereigns had new dams and aqueducts constructed and cared not to deprive the city of water. With the fountains, water dispensers, Maksem, cisterns and aqueducts constructed during Ottoman period, water gained its freedom and İstanbul met with water, water met with İstanbul. Today, İSKİ continues to provide safe, clean, and healthy water of good quality to İstanbul which has been the capital of water civilization.

MESSAGE FROM ISKI

Water has an important value in our lives throughout the history. Being well aware of this importance, previous civilizations and our civilization have been built on water. Cisterns, water arches and canals are the main indicators of water's importance. In our city, Istanbul, İSKİ keeps this tradition alive and lives through it. With traditional experience, today's technology, future's investments and personnel, İSKİ is the guaranteed provider of clean and healthy water.

By using its own technology and laboratories, with its professional personnel, İSKİ is building not just today's Istanbul, but also future's Istanbul for next generations. While Istanbul is getting ready for European Capital of Culture - 2010, İSKİ is endeavoring to complete its job in this manner. For this reason, İSKİ started working to spread quality system in 2005. Teaming up with Marmara Research Center (TÜBİTAK), universities and other related institutes, İSKİ is working on project progress and infrastructural systems.

As result, Water Quality Control Department is serving the customers with EN ISO/IEC 17025 (General requirements for competence of test and calibration laboratories) standards.

İSKİ is playing a key role in developing water quality control laboratories based on national sustainable development aims and technological progress in the world. In this manner İSKİ is helping country's economy by researching and applying projects which help country's short and long term needs and Turkey's integration with European Union.

We are pleased to share the information about our Water Quality Control Laboratory. As you will see, the water goes through many important and serious treatments before it reaches your tap. I wish your confidence in water will grow as it was before.

Best Regards

Mevlüt VURAL
İSKİ General Manager

CONTENTS

CHAPTER 1

THE ADDRESS OF WATER AND QUALITY: İSKİ

CHAPTER 2

THE ADVENTURE OF WATER POURING FROM FOUNTAINS

CHAPTER 3

PROPERTIES OF DRINKING WATER

CHAPTER 4

LABORATORY STUDIES

STUDIES OF THE BACTERIOLOGY LABORATORY

STUDIES OF THE CHEMISTRY LABORATORY

CHAPTER 5

WATER QUALITY CONTROL LABORATORY HAS BEEN ACCREDITED

OUR QUALITY VISION

WATER QUALITY LABORATORY

İSTANBUL BÜYÜKŞEHİR BELEDİYESİ

İ S K İ
İSTANBUL SU VE KANALİZASYON İDARESİ